

OPHIDIA II

Heart Stabilization System




Made in Japan

VITAL

OPHIDIA II

Stabilizer



1 Easier handle close

Developed design enables the arm can be fixed with less torque.

2 Arm tension indicator

Indicator displays arm tension, allowing for precise and consistent adjustment.

3 Fits various retractors

Adaptable to a range of retractors, from those used for midline sternotomy to small incisions.

4 Increased flexibility and better stabilization

The internal structure has been optimized to strengthen the arm and stabilize the anastomosis site, while enhancing flexibility for easier guidance of the foot step.

5 Tubing holder


Tubing holder maintains an organized working space.

6 Malleable pods for maximum anastomosis stability

The product design has been enhanced to improve heart surface adhesion.

OPHIDIA II

Positioner



7 Suspension maintains stable hemodynamics

Suspension absorbs and moves in sync with the heartbeats, allowing precise heart positioning while minimizing tissue trauma.

8 Reinforced and flexible arm

1 cm longer arm provides better visualization for heart rotation.

9 Suction cup with advanced vacuum stabilization

The tissue-conforming suction cup design ensures a secure lift and hold of the heart.


10 OPHIDIA suture holder

Suture holder designed specifically for OPHIDIA series securely places tension on the sutures.

※For more information on the types of retractors compatible with this suture holder, please feel free to contact us.

OPHIDIA II

OPHIDIA suture holder




OPHIDIA II tubing set



※Yellow tube is for a stabilizer and the blue one is for a positioner.

OPHIDIA II

Stabilizer SHORT SIZE / ORBIS



※All disposable items

Catalog #	Description	Units/Box
VTS-101R	OPHIDIA Stabilizer (Regular size)	1
VTS-102S	OPHIDIA Stabilizer (Short size)	
VTS-104R	OPHIDIA Stabilier ORBIS	
VTP-101R	OPHIDIA Positioner	6
VTB-001R	OPHIDIA Suture holder	



VITAL Corporation

Riverside Shinagawa Konan Bldg. 6F
1-2-5 Higashishinagawa, Shinagawa-ku,
Tokyo, 140-0002, JAPAN
TEL: +81 3 3458 1265 FAX: +81 3 3458 1263
E-mail: marketing@vital-j.co.jp

